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Founded in 2005, Harbin Hapro Electric Limited is a high-tech company located in National Universities Science Park in Harbin University of Science and Technology, Harbin, China. The products from Hapro are instruments, equipment which related to polymer industry and cable industry. Hapro keeps a good relation with Harbin University of Science and Technology. Joint working together with professors in the university, Hapro can also provide project solutions in cable industry and polymer industry. The products from Hapro Electric include Torque Rheometer System. Polymer Film Quality Inspection System. UV-Series LED Irradiation Cross-linking Equipment, etc. The products from Hapro were not only sold domestic (Hong Kong and Taiwan) but also distributed all over the world, Russia, Belarus, Ukraine, the United States, Germany, Canada, Australia, Italy, Hungary, South Africa, Panama, Mexico, Argentina, Brazil, Chile, Colombia, Egypt, Saudi Arabia, Turkey, the Philippines, India, Iran, Malaysia, Thailand, Singapore, South Korea, Japan, Vietnam and other more than 20 countries .

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Since the foundation of Harbin Hapro Electric Technology Co., Ltd. We have been doing our best to develop technology and products. Hapro is willing to provide high quality products by the professional and technical ability and bring more benefits by the high-tech productions to customers. Hapro looks forward to cooperating with you to build a brighter future together.

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RM-200C Torque Rheometer System consist of a main unit and function units, generally assembled by block models. The function units include the Mixer, Single Screw Extrusion unit, Conical Twin-screw Extrusion unit, Parallel Twin-screw (co and counter rotating) Extrusion unit, Rubber Extrusion unit. Downstream auxiliary equipment include Laboratory Film Blowing Machine, Laboratory Film Casting Machine, Laboratory Film Calender Machine, Film Quality Inspection System, High-pressure capillary rheometer, Wire coating machine, Two-Roller Mills, Platen Presses, Pelletizer, Dynamic testing oven, Sample cutter and Belt conveyor, etc. System can not only simulate actual process of polymer, evaluate machining performance, optimize process parameters, but also completing simulation experiment to match different die heads and auxiliary equipment.

RM-200 C 主机平台技术参数 / Main unit technology parameters			
电机功率 Power of motor	5.0/7.0 kW	动力输出转速 Speed range	0.1~150 / 200 / 250 r/min
转矩测量量程 Torque range	0 ~ 200 / 240 / 320 Nm	转速控制精度 Speed accuracy	0.05%F.S
转矩测量精度 Torque accuracy	0.1%F.S	动态温度控制精度 Dynamic controlling accuracy of temp	± 0.5° C
转矩测量方法 Torque measuring method	悬臂式动态扭矩测量 Cantilever	熔体压力测量量程 Melt pressure range	0.1 ~ 100 Mpa (规格可定制) specifications customizable

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The Hapro Mixer of Torque Rheometer system serve for testing the rheological behavior properties and extrusion ability of polymers and for analyzing and emulating production process in the laboratory. It is an ideal instrument and system in research and development as well as in practical application, for searching the material liquidity, plasticizing, heating and shear stability, etc. Save the time of practical production and the quality of material.

技术参数 Technology parameters				
	Mix-20	Mix-60	Mix-200	Mix-300
容积 Chamber volume	20ml	60ml	200ml	300ml
转子类型 Kinds of rotor	Roller 标准 Standard (Banbury、Cam、Delta、Sigma 可选 Available)			
转子速比 Speed ratio	3: 2			
最大扭矩 Max.torque	80Nm	200 / 240 / 320 Nm		
温控系统 Temp control system	RKC system			
温度传感器 Temp sensor	J 型热电偶 (K 型及其它可选) / J type thermocouple			
加热方式 Heating	3 区，电加热 / 3 Zones, Electric heating			
冷却方式 Cooling	空气冷却 Compress air cooling			
最高温度 Max.temp	350℃			
最大转速 Max. rotor speed	100 r/min	150 or 200 or 250 r/min		
压料方式 Loading	手动 (气动可选) Manually+Pneumatic			
安全设置 Security settings	电气 + 机械，双重保护 Electricity+Machinery+Safety key			

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Mixer of Torque Rheometer System is constituted by main unit and mixer unit. The mixer chamber is consist of three sections which are independently temperature controlled. There is a pair of rotors in it. Counter-rotating rotors provide high strength mixing and shearing action. Customer can choose the different chamber volumes, the different rotor types and loading ways to implement different applications. For teaching and demonstration, the visualization window can be configured.

● 转子类型
Kinds of Rotors

图例 Photo	名称 Name	功能 Function
	Roller	用于混合热塑性材料，如聚烯烃、聚氯乙烯、工程塑料等。 It is usually used in mixing thermoplastic materials, for example: PO, PVC, Engineering plastics, etc.
	Banbury	用于混合热塑弹性体，橡胶行业中应用更为普遍。 It is usually used in mixing thermoplastic elastomer and rubber material.
	Delta	用于热固性材料的混合和交联。 It is usually used in mixing and cross-linking thermoset polymer.
	Cam	用于混炼热塑料、陶瓷复合物、食品等粘稠高扭矩的物料。 It is usually used in mixing thermoplastic polymer, ceramic composites, food and other high-torque materials.
	Sigma	用于食品、药品和塑料溶胶等扭矩较小物料的混合。 It is usually used in mixing food, medicine and plastisol materials which have low-torque.

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
可视化视窗
Visual window for mixing


制样机
Sample tablet tool for
PVC thermal stability test


- 质量控制和原材料检测;
Quality control and testing of raw materials;
- 产品开发和配方优化;
Product development and recipe optimization;
- 加工性能研究与工艺优化;
Investigation and optimization of processing conditions;
- 共混料的流变特性测试;
Measurement of the rheological properties of blends and compounds;
- 热固性材料的交联(固化)特性测试;
Measurement of cross-linking and curing process of thermoset resin;
- 教学科研应用。
Applications in the field of teaching and research.

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
Mixer Test Report

熔合实验
Fusion test report

交联实验
Cross-linking test report

叠加实验
Compare test report

热稳定实验
Thermal stability test report

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The Conical Twin-screw Extrusion unit

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Laboratory Conical Twin-screw Extruder can be integrated by main unit and conical twin-screw extrusion unit. That can be used to simulate real production process in the laboratory. Because of its special structure, the residual material can be cleaned quickly. In the measuring extrusion unit, the sample material begin plastify under practice-oriented conditions and extruded through the die head. Matched well with different die heads, it's used for extruding the Belt, Tube, PCPO and granulation.

技术参数 Technology parameters			
螺筒分离结构 Barrel structure	拉出回转式 Pull and turn	螺杆类型 Screw type	普通型或双锥型 Normal or Biconical shape
螺杆材质 Screw material	38CrMoAlA 渗氮 or SKD11 渗钛	螺杆直径 Screw diameter	20-40mm
螺杆长度 Screw length	325mm	螺杆转速 Screw speed	200r/min
加热方式 Heating	电加热 5 区 双 PID Electrical heating,5 zones, Double PID	温控系统 Temp control system	RKC system
温度传感器 Temperature sensor	J 型热电偶 (K 型及其他可选) J type thermocouple	冷却方式 Cooling	压缩空气 Air compression
加热总功率 Heating power	3.5kW	压力范围 Pressure range	0-100Mpa (30Mpa 可选)
最高温度 Max. temp	350℃	产量 Output	1-10kg/hr

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Extrusion unit structure features


- 拉出回转式，便于清洁、分离与组装；
It's convenient to separate and assemble;
- 高精度滑轨以及旋转支撑，设备高稳定性；
High precision slides and bearing ,the higher equipment stability;
- 清洁时电气不分离，加热不间断；
No electric separation and heating uninterrupted;
- 防尘密封性好；
Dust proof sealed well;
- 清料托盘架自动弹出。
Automatic pop-up the tray rack.

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Application of The Conical Twin-screw Extrusion unit

- 质量控制和原材料检测；
Quality control and testing of raw materials;
- 配方开发与优化；
Testing the processing behavior during recipe development and optimization;
- 不同物料的挤出加工性能测试；
Testing of the extrusion ability of different materials;
- 加工工艺参数优化；
Investigation and optimization of processing conditions;
- 混料的流变性能和黏度性能测量；
Measurement of the rheological or viscosity properties of blends and compounds;
- 挤压降解产品；
Extrusion of degradable products;
- 小规模生产。
Mini-test/production.

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Laboratory Single Screw Extruder can be integrated by main unit and single screw extrusion unit. The screw with different structure,diameter, length-diameter ratio and compress ratio can be chosen depending on the sample material. Integrating different extruding dies and other auxiliary equipment, the system could simulate many functional experiments in the laboratory.

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Equipped with capillary dies, the extrusion unit could be used to measure the melt viscosity over a usual shear rate range. This so called extruding capillary Rheometer could obtain the viscosity information closed to extruding process. The melt flow and die pressure can be obtained under different rotation speed, by changing the screw speed and die cores with different length-diameter ratio, the relation between shear rate and melt viscosity can be calculated with software after dies entrance pressure calibration and nonlinear calibration.

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Laboratory Rubber Extruder can be integrated by main unit and rubber extrusion unit. Mainly used for the extrusion of rubber, thermosetting plastic and chemical raw materials.Integrating different extruding screw and dies heads, the system could simulate rubber products production process in the laboratory.

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The Parallel CO-rotating Twin-screw Extrusion unit

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Parallel twin screw extrusion unit (co-rotating), widely used in rubber and engineering resin blends, dispersion, modification, enhancement, chlorinated polypropylene and high water absorbing resin devolatilization, dedicated to the extrusion forming process, for example the blending and granulation experiments.

技术参数 Technical parameters

螺杆直径 Screw diameter	20 mm	长径比 Screw L:D	40:1/48:1 (其他规格可定制) specifications customizable
加热方式 Heating	电加热, 8区, 双PID Electrical heating, 8 Zones, Double PID	排气 Exhaust	真空排气 Vacuum
温控系统 Temp control system	RKC system	加热总功率 Heating power	6.4 kW
最高温 Max. temperature	450°C	螺杆转速 Max. speed of screw	480 r/min
温度传感器 Temperature sensor	J型热电偶 (其它型号可选) J type thermocouple	冷却方式 Cooling	压缩空气 Air compression
产量 Output	1 - 15 kg/ hr	压力范围 Pressure range	0-100 Mpa (30Mpa、50MPa、70MPa 可选)

同向平台结构特点 Structure features of co-rotating extrusion unit

- 螺杆元件采用“积木式”结构, 可根据物料和加工工艺的不同, 任意互换、组合使用;
Screw element use the "building block" structure, interchangeable and combination according to the different materials and processing technology;
- 筒体选用氮化钢、双金属材料, 耐磨、耐腐蚀, 延长使用寿命;
Barrels choose bimetallic material, wear-resisting, corrosion resistance and increasing the service life;
- 物料分散均匀, 混炼塑化效果好, 物料滞留时间短, 输送效率高;
Evenly material dispersion, better mixing plasticizing effect, short material residence time, high transmission efficiency;
- 采用加热冷却双PID, RKC温控系统, 反应灵敏, 温控精度高;
RKC temperature control system has a sensitive response and high precision of temperature control;
- 喂料方式多种可选, 可适用于颗粒料、粉料、液体料。
A variety of optional feeding way according to the particles, powder, liquid material.

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Loss-in-weight feeder system is used for the compounds of polymer materials and auxiliary formulations by medium and small parallel co-rotating twin-screw extruder in laboratory. The system is composed of several small-flow loss in weight feeder. The range of the weighing sensors and the type of feeding screws in the weight sensor scales can be combined according to the actual needs of the feeding, to achieve high-precision flow control of different types of materials and a variety of flow ranges, and the flow of all weight loss scales can be adjusted in tandem to achieve real-time setting of the mixing ratio of different compounds of the material formula.

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Melt pump is mainly used in high temperature and high viscosity polymer melt transportation, pressurization, metering. The melt pump pressurizes and stabilizes pressure of hot melt polymers coming from the extruder, then send them into extruder dies.

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The Parallel counter rotating
Twin-screw Extrusion unit

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The parallel twin screw extrusion unit (counter rotating) has a better plasticizing effect compared to the conical twin screw extruder. It is widely used in forming process of thermoplastic materials, for example the extrusion of sheet and pipe, or granulation for PVC resin.

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


锥双挤出报告
Conical twin-screw extrusion unit test report


平双挤出报告
Parallel twin-screw extrusion unit test report

材料剪切速率与粘度关系曲线
Shearing rate and viscosity curve

两种材料的粘度曲线对比
Compare between two different material' s viscosity curve

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Software of The Torque
Rheometer System

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

Torque rehometer simulation system

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Lab formulation blending extrusion products forming and on line examination system: The Hapro Torque Rheometer system, consisting of a main control unit and all kinds of function units. These units can be assembled free to constitute some fully functional simulation system of the practical production process, according to the different production and technique demands of polymer. For example, formulation, blending, extrusion, products molding, on line examination, etc.

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The Hapro simulation system, consisting of a main control unit, parallel twin-screw extrusion unit and melt pump based on formulation blending and extrusion with volume feeding, final products molding machine, such as Film casting machine, Film calendering machine, Film blowing machine, Wire coating machine, Pelletizer, Sample cutter, etc.

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Laboratory Film Calendering system adopts precision extrusion and high-performance tape calendering technology. The polymer material is extruded by conical-twin extrusion unit, through into heating calender rollers to get film. The equipment use oil heating method by mould temperature controller. Calender system is suitable for polymer materials film, such as PE, PP, PS, ABS, PVC, PET, PC, etc. A film quality inspection system can be equipped on the traction and take-up unit, used for inspecting the quality defects.


Table with 4 columns: Component, Specification, Component, Specification. Rows include Rolling structure, Diameter of roller, Roller width, and Distance between rollers.

Table with 4 columns: Component, Specification, Component, Specification. Rows include Heating method, Max. temp of roller, Power of mould temp controller, and Tubing connection.


Table with 4 columns: Component, Specification, Component, Specification. Rows include Pull and charge motor, AC motor, Charge control, and Tension adjustment.

触摸屏操作界面
Human machine interface

- 可设置手动模式与联动模式;
Manual mode and linkage mode;
- 设置压辊之间的转速比引离辊与牵引辊之间膜速比;
Setting the speed ratio between the roller , speed radio between the doff roller and pull roller;
- 实时监控回油温度;
Real time monitoring the feed back oil temperature;
- 可扩展视觉测量单元对膜进行检测;
Extensible visual measurement unit to test the film quality.

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Laboratory Film Calendering system, consisting of a main control unit, a conical twin-screw extrusion unit or single screw extrusion unit , a film calendering machine, a traction and take-up unit, on line inspect and detect unit. Profit from the precision extrusion and high-performance tape calendering technology. This system can be simulation of the practical production process.

应用领域
Application Field

可用于演示压延的工艺过程;
Demonstrating technological process of calendering;

材料配方研究与发展;
Material formulation research and development;

模拟压延实际生产过程;
Production process simulation;

聚合物材料质量控制。
Quality control of compounds and polymers.

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Film casting system adopts precision extrusion and high-performance tape casting technology. The extrudate leaving the Film casting die head can be made a thin film or sheet through forming and cooling. Suitable for polymer materials such as PE, PP, PS, ABS, PVC, PET, PC, etc.

技术参数 Technical parameters	
流延膜模具 Film casting die	150 × 1 mm, 1.8 kW (0-2.5mm 厚度可调 adjustment)
成膜厚度 Film thickness	0.02-0.5mm
成膜宽度 Max. film width	120mm
流延辊直径 Diameter of upside roller	200mm
流延调速 Roller speed	0-110mm/s
卷取方式 Crimp ways	张力控制 Tension adjustment
冷却方式 Cooling	水冷 Water
自动切边装置 AUTO cutting-edge device	

实验室塑料吹膜机
Laboratory Film Blowing Machine

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Thermoplastics Polymer have been melt and quantitative extrusion by the film blowing die in the extrusion unit. Then through ring cooling wind, blowing the herringbone plate, traction roller, coiling roller. Eventually made into film.

技术参数 Technical parameters	
吹膜模具 Film blowing die	Φ25 / Φ30 / Φ40 mm (规格可定制 specifications customizable)
模口间隙 Gap of die lip	0.6 / 1.0 mm
吹膜调速 Roller speed	0-140 mm/s
成膜折径 Film fold diameter	50-120 mm (规格可定制 specifications customizable)
成膜厚度 Film thickness	0.03-0.3 mm
卷取方式 Take off	重力收取 Gravity charge

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


The measurement in the shear rate region over app 10^{-2} – 10^6S^{-1} , which allows the simultaneous measurement at three different die lengths. The multi-functional equipment generation with unique characteristic data and highest possible precision, applied for quality control or research and development – a complete program range, covering nearly all needs!

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
Mainly used for testing the thermal stability of polymer powder, pellet, tubular, laminated samples. Dynamic testing oven can send the sample continuous or discontinuous, Experimenter could observe the process of sample color changes under different aging time. Testing oven can save the two thirds of test time compared to the traditional oven and turn the sample length into the time scale. Evaluating thermal stability properties of materials by contrasting the color change with color meter, Further improve the production process and optimize the product formula, etc.

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When an industrial high-speed and high-performance electronic camera was equipped on the film casting unit or a chill roller unit, a film or tape quality inspection system was constructed. This system can be used to inspect the foreign contamination particles and “fish eye” in the materials. This inspection is important to some high quality membrane material, for instance the high voltage cable insulation materials.

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压延式杂质检测仪
Contamination Detector For Film Calendaring

流延式杂质检测仪
Contamination Detector for Film Casting

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Characteristics of Contamination Detector

- 高性能，60MHz输出速度，线扫描工业相机；
High performance, up to 60MHz output, line scan industrial camera;
- 高速PCI数据采集系统；
High speed PCI data capture system;
- 亮度可调高亮LED线光源；
High light density LED linear light source(density adjustable);
- 张力反馈控制收取电机。
Tension feedback used to control crimp.

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Surface finish detection system of high-voltage
cable semiconductive material strip

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The structure of XLPE high voltage cable insulation is composed of an insulative layer and semiconductive layers. The semiconductiv material is a blend of polyolefin elastomer material and carbon black, which is used for the boundary layer of the electric field of the cable insulation layer. Since the physical properties of the semiconductive layer and the insulative layer are similar, they are also co-extruded in the production process, there is no air gap at the interface. The electric field in the insulation from the high potential of the core shielding layer to the low potential of the external shielding layer is perfect in XLPE insulation layer, which can withstand higher operating voltage and service life. According to the principle of higher electric field strength at the small curvature radius of the equipotential surface, the semi-conductive layer faces the insulation. If there are protrusions on the interface, it will cause the increase of the electric field strength to damage the insulation. Therefore for, the production of semi-conductive materials and cable insulation systems, it is required to limit the number and size of the protrusions on the surface of the semi-conductive layer strictly. The specific implementation of the quality control system is to control the number and size of the protrusions on the surface of the semiconductive material strip.

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The development of high-voltage cable semi-conductive shielding material surface finish detection system is based on low-temperature extrusion and casting technology of polymer materials, and electronic camera optical detection technology (or machine vision technology). The system is equipped with a high-speed vision camera and a special LED line-light source to detect size information of the protrusion in real time. The image processing capacity is 15 million pixels per second. It can achieve high-speed and high-resolution measurement. At the same time, the camera system has a wide field of view and high pixels. It can synchronous measure the calibration area while measuring the sample to calibrate real-time.

挤出系统 / Extrusion system	
电源 Power supply	3×380V,50/60 Hz ,+N +PE
电机 Motor	施耐德伺服电机 Schneider servo motor
电机功率 Power of motor	1.8 kW
减速比 Reduction ratio	1:35
螺杆直径 Screw diameter	20mm
螺杆长径比 Screw L:D	25 :1
螺杆材质 Screw material	进口高强度合金钢 Imported high-strength alloy steel
螺杆压缩比 Screw compression	1.5:1
温控系统 Temperature control system	日本 RKC, 加热冷却双 PID 控制 RKC company of Japan, heating and cooling dual PID control
喂料方式 Feeding	体积喂料 Volumetric feeder
最高可控温度 Maximum controllable temperature	350℃
温度传感器 Temperature sensor	J 型热电偶 J type thermocouple
加热区 Heating zone	4 区加热 (螺筒 3 区加热冷却, 模具 1 区加热) 4 zone heating (3 zone heating and cooling of the screw barrel, heating zone 1 of the die)
加热方式 Heating	电加热 Electrical heating
加热总功率 Heating power	4.0kW
冷却方式 Cooling	风冷 Wind
模具 Die	口模尺寸 60mm×1.5mm lip of die 60mm×1.5mm

成膜系统 / Film Casting system	
电源 Power supply	3×380V,50/60 Hz ,+N +PE
流延辊尺寸 Size of casting roll	Φ200mm×200mm
流延辊中心高 Center height of casting roller	1022 mm
流延成膜宽度 Reduction ratio	25~30 mm
流延成膜厚度 Thickness film of casting	0.5mm±0.02
流延辊控温方式 Method of casting roll temperature control	冷热一体机恒温控制 constant temperature control of the cooling and heating integrated machine
流延辊驱动电机 Casting roller drive motor 功率 power	伺服电机 servo motor 640W
牵引辊驱动电机 Traction roller drive motor 功率 power	伺服电机 servo motor 640W

光学检测系统 / Optical inspection system	
测量厚度 Thickness of measurement	0.4~0.5 mm
凸起高度分辨率 Resolution of protrusion height	20 μ m
视场宽度 Width of field	0~75 mm
检测宽度 Width of detection	44mm
最小成像宽度 Width of the minimum imaging	10 μ m
最小测量尺寸宽度 Width of the minimum measurement size	30 μ m
检出率 Detection rates	100%
测试速度 Speed of measurement	>1kg/h
室内环境 Work environment	一万级 The class of ten thousand
使用温度 Operating temperature	0 ~ +40℃
收膜速度 Film-drawing speed	0~110mm/s
环境湿度 Environment humidity	< 65%
测量辊参数 Measuring roller parameters	
直径 Diameter	40mm
同轴精度 Coaxial accuracy	3 μ m
测试光源 Test light source	220mm 长线扫描光源 220mm line-scanning light source

绕包模塑式高压电缆软接头
绝缘带和屏蔽带制备系统
Insulating tape and shielding tape production
system for molded high-voltage cable flexible joint

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The preparation method of the wrap-around molded flexible joint is complementary to the extrusion molding method. It is very effective for the preparation of the high-voltage cable flexible joint. Especially when the fault is repaired, the wrap-around molding method is more convenient and effective. In addition, foreign reports indicate that all cable joints below 138kV tend to use wrap-around molding technology, which is more meaningful for domestic system applications below 110kV. The insulation and shielding layer wrapping tape production system are developed based on the integrated technology of polymer material extrusion, casting and slitting, equipped with a real-time detection system for impurities in the insulation tape with electronic cameras.

Insulation tape production system composition: defect-free single screw extruder, casting and slitting machine, film casting die, electronic camera impurity inspection, dual-channel independent mold temperature machine.

挤出系统 / Extrusion system	
电源 Power supply	3 × 380V,50/60 Hz ,+N +PE
电机 Motor	3 kW 伺服电机 3kW servo motor
减速比 Reduction ratio	1:15
螺杆直径 Screw diameter	30mm
长径比（L:D） Screw L:D	25 :1
螺杆材质 Screw material	38 CrMoAlA
温控系统 Temp control system	日本 RKC, 加热冷却双 PID 控制 RKC company of Japan, heating and cooling dual PID control
最高可控温度 Maximum controllable temperature	350℃
温度传感器 Temperature sensor	J 型热电偶 J type thermocouple
喂料方式 Feeding	体积喂料 Volumetric feeder
加热区 Heating zone	8 区电加热（螺筒 3 区加热冷却，机颈 1 区，模具共 4 区加热） 8 zones of electric heating (3 zones for heating and cooling of the screw barrel, 1 zone for the machine neck, 4 zones for the die heating)
加热总功率 Heating power	6.0kW
冷却方式 Cooling	风冷 Wind

流延模具 / Film casting die			
模具总体尺寸 Overall die size	320 × 250 × 170mm	模口宽度 Width of die lip	250mm
模具模口可调尺寸 Adjustable size of die lip	0.02~2.5mm	成膜厚度 Thickness of film	0.02~0.5mm

流延分切系统 / Casting slitting system			
电源 Power supply	3 × 380V,50/60 Hz ,+N +PE	流延辊直径 Casting roller diameter	150mm
流延辊中心高 Center height of casting roller	1022 mm	流延膜模具尺寸 Film casting	150mm*0.02~2.5mm 厚度可调 Adjustable thickness
流延成膜宽度 Film width of casting	< 120 mm	流延成膜厚度 Film thickness of casting	0.08~0.5mm
流延辊冷却方式 Casting roller cooling method		水冷（水温可控） Water(controllable water temperature)	

杂质检测光学系统 / Impurity detection optical system			
测量厚度 Measurement thickness	0.04~1 mm	分辨率 Resolution	20 μ m
视场宽度 Field width	0~100 mm	检出率 Detection rates	100%
测试速度 Speed of measurement	>1kg/h	室内环境 Work environment	一万级 The class of ten thousand
使用温度 Operating temperature	0 ~ +40℃	收膜速度 Charging film speed	0~110mm/s
环境湿度 Environment humidity		< 65%	

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
The system consists of a single screw extruder unit, a melt pump and a pressure filter die. The extrusion principle is a test mixture composed of a color masterbatch and a basic sample polymer. Through an extruder equip a melt pump and a plate filter, a melt pressure sensor is installed in front of the filter, and the filter value is calculated by the difference between the initial pressure and the maximum pressure to evaluate the dispersion performance of the colorant. In order to ensure that the melt pump is completely full and the melting is completely uniform, the system has the function of pressure feedback.

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Experimental simulation line of high temperature 3D printing wire extrusion

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The system consists of a high temperature extrusion unit, a cooling setting unit, a diameter measuring unit and a traction collection unit. The maximum extrusion temperature of the high temperature extrusion unit is 400 °C. The extrusion temperature and screw speed are controlled by the touch screen, and the changes of temperature, speed and pressure in the extrusion process can be displayed in real time. Taking the measured real-time diameter as the target, the traction speed is controlled by the diameter measuring device to extrude the wire with the diameter of 1.75mm or 3mm, and the diameter error is less than 0.05mm.

高温单螺杆挤出机 / High temperature single screw extruder			
驱动方式 Driving mode	伺服驱动 servo drive	电机功率 Power of motor	1.5kW
减速器减速比 Reducer deceleration ratio	15:1	螺杆转速 Screw speed	95r/min
加热区 heating zone	5 区 5 zone	温度范围 Temperature range	室温 ~400℃ room temperature ~400℃
温控精度 Temperature control precision	±0.5℃	压力测量范围 Pressure measuring range	0.1~100Mpa
螺杆直径 Screw diameter	20mm	长径比 Screw L:D	30:1
压缩比 Compression ratio	3:1	冷却方式 Cooling	喂料口水冷，螺筒气冷 water cooling for feed zone, air cooling for barrel

线材模具 / Wire die			
线材模具 Wire die	Φ2.25mm	模具处理 Die treatment	内表面镀铬并镜面抛光 hard chrome plating on the inner surface and mirror polishing
制品规格 Product specification	直径 DN1.75mm diameter DN1.75mm	制品公差 Product tolerance	±0.05mm

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Planetary mixer is suitable for the reaction, mixing, dispersing and blending of materials such as liquid and liquid, solid and liquid, paste and paste, powder and paste, as well as materials with a wide disparity in density, from low viscosity to high viscosity. It is composed of planetary mixing blade and scraper.

Features: Full mixing, simple structure, easy disassembly and assembly, different mixing purposes can be met by replacing the paddle type, and the effect is better than ordinary mixers.

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
The vacuum pipe (Micropipe) forming system is composed of torque rheometer main unit, single screw extrusion unit, vacuum pipe forming unit, diameter measuring unit and Belt Conveyor. Extruder can be equipped with high temperature and corrosion resistant materials, suitable for PP, PE, PVDF and other materials .It can simulate pipe extrusion forming intuitively and effectively.

可移动式测控主机平台 / Movable Platform of Main unit	
性能指标 Performance Index	性能参数 Performance Parameter
电机功率 Power of Motor	5.0 kW
电机速度 Motor Speed	3000 r/min
减速比 Reducer ratio	1:15
动力输出转速 PTO speed	0.1~200 r/min
转速精度 Speed Accuracy	0.05% F.S
扭矩范围 Torque Range	0~240 Nm
扭矩精度 Torque Accuracy	0.1% F.S
静态温度控制精度 Static Measurement Accuracy of Temp	0.1℃
动态温度控制精度 Dynamic Measurement Accuracy of Temp	+0.5℃
熔体压力量程 Melt pressure Range	0.1~100 Mpa
压力精度 Melt pressure Accuracy	0.5% F.S



可移动式挤塑平台 / Movable Platform of Single screw Extruder	
性能指标 Performance Index	性能参数 Performance Parameter
螺杆结构 Screw configuration	等距变深 Isometric darker
螺杆长径比 Screw L:D	30:1
螺杆直径 Screw Diameter	20 mm
螺杆压缩比 Compression Ratio	3.6:1
温控系统 Temp control system	日本 RKC / RKC from Japan
加热 / 螺筒 5 区 Heating	电加热 Electrical Heating / barrel 5 Zones
加热功率 Heating power	4.6 kW
最高温度 Max.Temperature	300℃
喂料方式 Feeding	体积喂料 Volumetric feeder

真空管成型设备 / Vacuum pipe forming unit	
性能指标 Performance Index	性能参数 Performance Parameter
水槽总长 Water trough length	1500mm
真空段长 Vacuum section length	500mm
水槽宽度 Water trough width	200mm
真空泵 Vacuum pump	水环式 Water ring
成管直径 Final pipe diameter	1.3~1.5mm
壁厚 Wall thickness	0.1mm

皮带式牵引机 / Belt Conveyor	
性能指标 Performance Index	性能参数 Performance Parameter
牵引有效长度 Effective hauling length	200mm
牵引有效宽度 Effective traction width	60mm
牵引速度 Traction speed	2~20m/min adjustment
收取方式 taken mode	自动收卷 Automatic winding
收卷线速度 Winding line speed	20m/min
适用管径规格 Applicable pipe diameter specifications	≤ Within Φ2 mm
测径仪测量范围 Measuring range of caliper	0~15mm
测量精度 Measurement accuracy	±0.03mm

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
Laboratory polymer extrusion molding test system

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Mainly used for extrusion testing of conventional and high filling UPVC formula sheet, pipe and profile. It can complete the simulation of extrusion production with conical twin screw.

可移动式主控平台 / Movable platform of Main unit	
性能指标 Performance Index	性能参数 Performance Parameter
电机功率 Power of Motor	11.0 kW
电机速度 Motor Speed	1500 r/min
减速比 Reducer Rate	1:10
速度范围 Speed Range	0.1~150 r/min
转速精度 Speed Accuracy	0.05% F.S
扭矩范围 Torque Range	0~900 Nm
扭矩精度 Torque Accuracy	0.1% F.S
静态温度控制精度 Static Measurement Accuracy of Temp	0.1℃
动态温度控制精度 Dynamic Measurement Accuracy of Temp	+0.5℃
熔体压力量程 Melt pressure Range	0.1~30 MPa
压力精度 Melt pressure Accuracy	0.5% F.S

可移动式锥形双螺杆挤出平台 / Movable platform of Conical Twin Screw Extruder	
性能指标 Performance Index	性能参数 Performance Parameter
螺杆组合方式 Barrel combination	气动分离螺杆筒，方便清理螺杆 Pneumatic pull and push barrel and convenient to clean the screw
双螺杆规格 Twin-screw specifications	前端直径 little-end:35 mm 后端直径 Big-end:81mm
螺棱长度 Screw arris length	865mm
温控系统 Temp control system	RKC from Japan
加热 Heating	Electrical Heating/8 Zones (Barrel 4,Die head 4)
螺杆加热功率 Barrel heating power	9 kW
模具 Die head	T- Die W×T : 100×4~6mm adjustable
最高温 Max.Temperature	350℃
冷却 Cooling	压缩空气 Air Compression
排气 Exhaust way	真空 Vacuum
压力量程 Pressure range	0~30 Mpa

板材成型实验机 PVC sheet forming machine	
性能指标 Performance Index	性能参数 Performance Parameter
辊轮规格 Roller specification	直径和长度 Φ110×190mm
辊转速 Roller rotation speed	0.1~40 rpm
辊温控 Roller temperature control mode	水温机 Mold temperature machine, medium is water
温度范围 Temp. Range	5~90℃
最大成型宽度 Max. Forming width	120 mm
成型厚度 Range of forming thickness	1~7.0 mm
牵引速度 Traction speed	0.1~10 m/min
冷却 Cooling	水冷 Water